

DNase Test Agar Plate w/ Methyl Green

MP1419

Intended Use:

Recommended for detection of deoxyribonuclease activity of bacteria and fungi and especially for identification of pathogenic Staphylococci.

Composition**

Ingredients	Gms / Litre
Tryptose	20.000
Deoxyribonucleic acid (DNA)	2.000
Sodium chloride	5.000
Methyl green	0.050
Agar	15.000
Final pH (at 25°C)	7.3±0.2

**Formula adjusted, standardized to suit performance parameters

Directions

Either streak, inoculate or surface spread the test inoculum (50-100 CFU) aseptically on the plate.

Principle And Interpretation

DNase test Agar is used for detecting deoxyribonuclease activity of bacteria and fungi and particularly for identification of pathogenic Staphylococci (1). DNase producing organisms exhibit clear zone around growth against green background. Reagent addition is not required (2). This medium is based on modification of the procedure for detecting DNase-producing bacteria as per Smith, Hanoch, and Rhoden (3) and Jefferies, Holtman and Guse (4). The medium supports growth of both gram positive and gram-negative bacteria.

Tryptose serves as nitrogenous source for the organisms. DNase produced by microorganisms depolymerizes the DNA substrate in the medium. Methyl green fades into a colourless compound producing distinct clear zones surrounding colonies (or band/spot inocula) in an otherwise green coloured medium. Methyl green requires a highly polymerized DNA substrate (5) and it combines with polymerized DNA forming a stable, green complex at pH 7.5 (6-8). As hydrolysis progresses, methyl green is released and when not combined at this pH it fades and becomes a colourless compound. Therefore clear zones are observed (7,9).

Type of specimen

Isolated Microorganism from non clinical samples

Specimen Collection and Handling:

For isolated microorganism samples, follow appropriate techniques for sample collection and processing as per guidelines (1,3). After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/ face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheet

Limitations :

1. The test organisms must be in pure culture and 18-24 hours old.
2. Individual organisms differ in their growth requirement and may show variable growth patterns on the medium.
3. Each lot of the medium has been tested for the organisms specified on the COA. It is recommended to users to validate the medium for any specific microorganism other than mentioned in the COA based on the user's unique requirement.
4. Further biochemical and serological test are need to be carried out for confirmation.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the expiry period when stored at recommended temperature.

Quality Control

Appearance

Sterile DNase Test Agar Plate w/ Methyl Green in 90 mm disposable plates with smooth surface and absence of black particles/cracks/bubbles

Colour of medium

Greenish blue coloured medium

pH

7.10-7.50

Sterility Check

Passes release criteria

Cultural Response

Cultural characteristics observed after an incubation at 35-37°C for 18-24 hours

Organism	Inoculum (CFU)	Growth	DNase Activity
<i>Serratia marcescens</i> ATCC 8100	50-100	luxuriant	positive, clear halo around the growth.
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> ATCC 25923 (00034*)	50-100	luxuriant	positive, clear halo around the growth.
<i>Staphylococcus epidermidis</i> ATCC 12228 (00036*)	50-100	luxuriant	negative reaction
<i>Streptococcus pyogenes</i> ATCC 19615	50-100	luxuriant	positive, clear halo around the growth.

Key : *Corresponding WDCM numbers.

Storage and Shelf Life

On receipt store between 2-8°C. Use before expiry date on the label. Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with sample must be decontaminated and disposed of in accordance with current laboratory techniques (10,11).

Reference

1. Macfaddin, J. F., 1985, Media for Isolation-Cultivation-Identification-Maintenance of Medical Bacteria, Volume1 Williams, Wilkins, Baltimore.
2. Schreier 1969. Am. J. Clin. Pathol. 51:711.
3. Jeffries C.D.; Holtman, D.F.; and Guse, D.G (1957) J. Bacteriol., 73, 590.
4. Smith, P.B., Hancock, G. A., and Rhoden, D. L (1969) Appl. Microbiol., 18,991.
5. Lachica, R.V.F. and Deibel, R. H (1969). Appl. Environ, Microbiol., 32 (4), 633.
6. Kurnick, N.B (1947). Cold Spring Harbor Symp. Quant. Biol., 12, 141.
7. Kurnick, N.B (1950) Arch. Biochem., 29, 41.
8. Kurnick, N.B and Foster, M. (1950) J. Gen. Microbiol. 33. 243.
9. Kurnick, N.B and Foster, M. (1950) J. Gen. Physiol. 34, 147.
10. Isenberg, H.D. Clinical Microbiology Procedures Handbook 2nd Edition.
11. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.

Revision : 01/2023

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.