

Perfringens Agar Plate

MP579

Intended Use:

Recommended for selective isolation and enumeration of *Clostridium perfringens* in food

Composition**

Ingredients	Gms / Litre
Tryptone	15.000
Soya peptone	5.000
Yeast extract	5.000
HL extract #	7.000
Ferric ammonium citrate	1.000
Sodium metabisulphite	1.000
Tris buffer	1.500
Agar	15.000
Final pH (at 25°C)	7.3±0.2

**Formula adjusted, standardized to suit performance parameters

Equivalent to Liver extract

Directions

Either streak, inoculate or surface spread the test inoculum (50-100 CFU) aseptically on the plate.

Principle And Interpretation

Clostridial species are one of the major causes of food poisoning/ gastrointestinal illnesses. They are gram-positive spore-forming rods that occur naturally in the soil (1). Foods commonly contaminated with *Clostridium perfringens* include meat, meat pies, poultry, stews and gravy. Among the family are: *Clostridium botulinum* which produces one of the most potent toxins in existence; *Clostridium tetani*, causative agent of tetanus; and *C. perfringens* commonly found in wound infections and diarrhoea cases. The use of toxins to damage the host is a method deployed by many bacterial pathogens. The major virulence factor of *C. perfringens* is the CPE enterotoxin, which is secreted upon invasion of the host gut, and contributes to food poisoning and other gastrointestinal illnesses (1).

Perfringens Agar (O.P.S.P.) is based on the formula developed by Handford (2) and is used as a selective medium for isolation and enumeration of *C. perfringens* in foods (3).

Tryptone, yeast extract, Soya peptone and HL extract supply most of the essential nitrogenous nutrients, vitamin B complex and trace ingredients for the growth of *C.perfringens*. Sodium metabisulphite and ferric ammonium citrate are used as indicators of sulphate reduction by *C. perfringens*, which produces black colonies. Tris buffer helps in maintaining buffering action. The antibiotics sulphadiazine, oleandomycin and polymyxin B make the medium highly selective inhibiting sulphite-reducing bacteria other than *C. perfringens* such as *Salmonella*, *Bacillus* species, *Proteus* species, *Staphylococci* etc.

Prepare 10 fold dilution of a 10 % homogenate of the food sample in 0.1 % Peptone Water (M028). Viable counts of *C.perfringens* bacilli or spores are obtained by plating 0.1 ml of different dilutions onto duplicate plates of blood agar containing 5 mg/lit of gentamicin/lit. Incubate at 37°C for 18-24 hours in two sets, one anaerobically and another aerobically. Alternatively incorporate 1 ml of the dilution into 25 ml of molten and cooled Perfringens Agar (O.P.S.P.) containing supplements. Incubate anaerobically for 18-24 hours at 37°C. Perfringens Agar with supplements gives high degree of selectivity and specificity.

Type of specimen

Food samples.

Specimen Collection and Handling:

Prepare 10 fold dilution of a 10 % homogenate of the food sample in 0.1 % Peptone Water (M028). Viable counts of *C. perfringens* bacilli or spores are obtained by plating 0.1 ml of different dilutions onto duplicate plates of blood agar containing 5 mg/lit of gentamicin/lit. Incubate at 37°C for 18-24 hours in two sets, one anaerobically and another aerobically. Alternatively incorporate 1 ml of the dilution into 25 ml of molten and cooled Perfringens Agar (O.P.S.P.) containing supplements. Incubate anaerobically for 18-24 hours at 37°C. Perfringens Agar with supplements gives high degree of selectivity and specificity.

Warning and Precautions :

Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/ face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheets

Limitations :

1. Further biochemical and serological tests must be carried out for further identification.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the expiry period when stored at recommended temperature.

Quality Control

Appearance

Sterile Perfringens Agar Plate in 90 mm disposable plates.

Colour of medium

Amber coloured medium

Quantity of medium

25 ml of medium in 90 mm disposable plates.

Reaction

7.10-7.50

Sterility Test

Passes release criteria

Cultural Response

Cultural characteristics observed after an incubation at 35-37°C for 18-48 hours

Organism	Inoculum (CFU)	Growth	Recovery	Colour of colony
<i>Bacillus subtilis subsp. spizizenii</i> ATCC 6633 (00003*)	$\geq 10^4$	inhibited	0%	
<i>Clostridium bifermentans</i> ATCC 17837	$\geq 10^4$	inhibited	0%	
<i>Clostridium butyricum</i> ATCC 13732	$\geq 10^4$	inhibited	0%	
<i>Clostridium perfringens</i> ATCC 12924	50-100	luxuriant	$\geq 50\%$	black
<i>Enterococcus faecalis</i> ATCC 29212 (00087*)	50-100	none-poor	$\leq 10\%$	white, if any
<i>Proteus vulgaris</i> ATCC 13315	$\geq 10^4$	inhibited	0%	
<i>Salmonella</i> Typhi ATCC 6539	$\geq 10^4$	inhibited	0%	
<i>Staphylococcus aureus subsp. aureus</i> ATCC 25923 (00034*)	$\geq 10^4$	inhibited	0%	

Key : *Corresponding WDCM numbers.

Storage and Shelf Life

On receipt store between 2-8°C Use before expiry date on the label.

Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with sample must be decontaminated and disposed of in accordance with current laboratory techniques (4,5).

Reference

1. Czeczulin J. R., Hanna P. C., McClane B. A., 1993, Infect. Immun. 61: 3429-3439.
2. Handford P. M., 1974, J. Appl. Bacteriol., 37: 559.
3. Hauschild A. H. W. et al, 1977, ICMSF Methods Studies VIII, Can. J. Microbiol., 23:884.
4. Isenberg, H.D. Clinical Microbiology Procedures Handbook 2nd Edition.
5. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.
6. Salfinger Y., and Tortorello M.L., 2015, Compendium of Methods for the Microbiological Examination of Foods, 5th Ed., American Public Health Association, Washington, D.C.

Revision : 00/ 2019

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.