

Plate Count HiVeg[®] Agar

MV091A

Intended Use:

Recommended for determining plate counts of microorganisms in milk and dairy products by pour plate technique.

Composition**

Ingredients	g / L
HiVeg [®] hydrolysate	5.000
Yeast extract	2.500
Dextrose (Glucose)	1.000
Agar	9.000
Final pH (at 25°C)	7.0±0.2

**Formula adjusted, standardized to suit performance parameters

Directions

Suspend 17.5 grams in 1000 ml purified/distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C. Mix well and pour into sterile Petri plates.

Principle And Interpretation

Plate Count Agar is equivalent to the medium recommended by APHA for the isolation of microorganisms in milk and other dairy products (1). Plate Count HiVeg[®] Agar is prepared by using vegetable peptones in place of animal based peptones which make the media free of BSE/TSE risks.

HiVeg[®] hydrolysate provides amino acids and other complex nitrogenous substances. Yeast extract supplies Vitamin B complex. APHA recommends pour plate technique. The samples are diluted and appropriate dilutions are placed in Petri plates. Sterile molten agar is added to these plates and plates are rotated gently to ensure uniform mixing of the sample with agar. Plate Count Agar is also used for the estimation of the number of live heterotrophic bacteria in water.

Type of specimen

Dairy samples

Specimen Collection and Handling:

For dairy samples, follow appropriate techniques for sample collection and processing as per guidelines (1,2). After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions :

Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/ face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheets.

Limitations :

1. Further biochemical and serological tests must be carried out for further identification.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the expiry period when stored at recommended temperature.

Quality Control

Appearance

Cream to yellow coloured homogeneous free flowing powder

Gelling

Firm, comparable with 0.9% Agar gel.

Colour and Clarity of prepared medium

Light yellow coloured clear to slightly opalescent gel forms in Petriplates.

Reaction

Reaction of 1.75% w/v aqueous solution at 25°C. pH : 7.0±0.2

pH

6.80-7.20

Cultural Response

Cultural characteristics after an incubation at 35 - 37°C for 18 - 24 hours.

Organism	Inoculum (CFU)	Growth
** <i>Bacillus spizizenii</i> ATCC 6633 (00003*)	50-100	luxuriant
<i>Escherichia coli</i> ATCC 25922 (00013*)	50-100	luxuriant
<i>Enterococcus faecalis</i> ATCC 29212 (00087*)	50-100	luxuriant
<i>Lactobacillus acidophilus</i> ATCC 4356 (00098*)	50-100	luxuriant
<i>Lactobacillus casei</i> ATCC 9595	50-100	luxuriant
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> ATCC 25923 (00034*)	50-100	luxuriant
<i>Streptococcus pyogenes</i> ATCC 19615	50-100	luxuriant

Key :* - corresponding WDCM numbers

Formerly known as *Bacillus subtilis* subsp. *spizizenii*.Storage and Shelf Life**

Store between 10-30°C in a tightly closed container and the prepared medium at 20-30°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition. Seal the container tightly after use. Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with sample must be decontaminated and disposed of in accordance with current laboratory techniques (3,4).

Reference

1. American Public Health Association, 1978, Standard Methods for the Examination of Dairy Products, 14th ed., APHA Inc. Washington, D.C.
2. Wehr H. M. and Frank J. H., 2004, Standard Methods for the Microbiological Examination of Dairy Products, 17th Ed., APHA Inc., Washington, D.C.
3. Isenberg, H.D. Clinical Microbiology Procedures Handbook 2nd Edition.
4. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.

Revision :01 / 2025

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.