

***HiCrome™ Enterobacter sakazakii Agar, Modified**
(HiCrome™ Cronobacter sakazakii HiVeg™ Agar, Modified)
(**Enterobacter sakazakii* now referred as *Cronobacter sakazakii*)

MV1641

Intended Use:

Recommended for the isolation and identification of **Cronobacter sakazakii* from milk and milk products.
Key: *: Formerly known as *Enterobacter sakazakii*

Composition**

Ingredients	g / L
HiVeg™ hydrolysate	7.000
Yeast extract	3.000
Sodium chloride	5.000
Synthetic detergent No. III	0.600
Chromogenic substrate	0.150
Crystal violet	0.002
Agar	15.000
Final pH (at 25°C)	7.0±0.2

**Formula adjusted, standardized to suit performance parameters

Directions

Suspend 30.75 grams in 1000 ml purified / distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C and pour into sterile Petri plates.

Principle And Interpretation

Enterobacter species are widely distributed in nature occurring in fresh water, soil, sewage, plants, vegetables, animal and human faeces. **Cronobacter sakazakii* has been closely associated with neonatal meningitis and sepsis (1). HiCrome™ Enterobacter sakazakii Agar, Modified is recommended by ISO Committee for the isolation and identification of *C. sakazakii* (2). HiCrome™ Enterobacter sakazakii HiVeg™ Agar, Modified is prepared by completely replacing animal based peptones with vegetable peptones to avoid BSE/TSE risks associated with animal peptones. The chromogenic substrate is cleaved specifically (3) by **C. sakazakii* resulting in the formation of blue green colonies. Other organisms, which do not cleave this substrate, produce colorless to slightly violet coloured colonies. HiVeg™ hydrolysate and yeast extract provide the essential growth nutrients along with nitrogenous and carbonaceous compounds. Sodium chloride helps in maintaining the osmotic equilibrium of the medium. Synthetic detergent No. III and crystal violet inhibits the accompanying gram-positive flora.

Key: *: Formerly known as *Enterobacter sakazakii*

Type of specimen

Food and Dairy samples; Water and wastewater samples

Specimen Collection and Handling

For food and dairy samples, follow appropriate techniques for sample collection and processing as per guidelines (1,4,5).
For water samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards (6).
After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions

Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheets.

Limitations

1. Some species may show poor growth due to nutritional variations.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the expiry period when stored at recommended temperature.

Quality Control

Appearance

Light yellow to pink homogeneous free flowing powder

Gelling

Firm, comparable with 1.5% Agar gel

Colour and Clarity of prepared medium

Light purple coloured, clear to slightly opalescent gel forms in Petri plates

Reaction

Reaction of 3.07% w/v aqueous solution at 25°C. pH : 7.0±0.2

pH

6.80-7.20

Cultural Response

Cultural characteristics observed after an incubation at 44±1°C for 18-24 hours.

Organism	Inoculum (CFU)	Growth	Recovery	Colour of Colony
# <i>Klebsiella aerogenes</i> ATCC 13048 (00175*)	50-100	good-luxuriant	≥50%	colourless with blue centre.
<i>Enterococcus faecalis</i> ATCC 29212 (00087*)	≥10 ⁴	inhibited	0%	
<i>Staphylococcus aureus</i> ATCC 25923(00034*)	≥10 ⁴	inhibited	0%	
<i>Escherichia coli</i> ATCC 25922 (00013*)	50-100	good-luxuriant	≥50%	colourless
## <i>Cronobacter sakazakii</i> ATCC 12868	50-100	good-luxuriant	≥50%	blue-green

Key : *Corresponding WDCM numbers. # - Formerly known as *Enterobacter aerogenes*

##-Formerly known as *Enterobacter sakazakii*

Storage and Shelf Life

Store between 15-25°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with sample must be decontaminated and disposed of in accordance with current laboratory techniques (7,8).

Reference

1. American Public Health Association, Standard Methods for the Examination of Dairy Products, 1978, 14th Ed., Washington D.C.
2. Muytjens H. L., Zanen H. C., Sonderkamp H. J. et al, J. Clin Microbiol 18:115-120, 1983.
3. International Organization for Standardization Draft ISO/ TS 22964, 2006 (E).
4. Salfinger Y., and Tortorello M.L. Fifth (Ed.), 2001, Compendium of Methods for the Microbiological Examination of Foods, 5th Ed., American Public Health Association, Washington, D.C.
5. Wehr H. M. and Frank J. H., 2004, Standard Methods for the Microbiological Examination of Dairy Products, 17th Ed., APHA Inc., Washington, D.C.

6. Lipps WC, Braun-Howland EB, Baxter TE, eds. Standard methods for the Examination of Water and Wastewater, 24th ed. Washington DC:APHA Press; 2023.
7. Isenberg, (2nd Ed.), Clinical Microbiology Procedures handbook, Volume 1, American Society for Microbiology, Washington, D.C.
8. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.

Revision : 05/2024

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.