

Yeast Malt HiVeg™ Broth (YM HiVeg™ Broth)

MV425

Intended Use:

Recommended for isolation and cultivation of yeasts, moulds and other aciduric microorganisms.

Composition**

Ingredients	g / L
HiVeg™ peptone	5.000
Yeast extract	3.000
Malt extract	3.000
Dextrose (Glucose)	10.000
Final pH (at 25°C)	6.2±0.2

**Formula adjusted, standardized to suit performance parameters

Directions

Suspend 10.5 grams in 490 ml purified/distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C. For preparing selective media acidify the media upto pH 3.0 to 4.0 by aseptically adding 1 vial of 10% Lactic Acid solution (FD095) or add antibiotics. **DO NOT HEAT** the media after addition of acid or antibiotics. Mix well and dispense as desired.

Principle And Interpretation

Yeast Malt Broth is formulated as per Wickerham (1,2) for isolation and cultivation of yeasts, moulds and other aciduric microorganisms. Fungistatic materials such as sodium propionate and diphenyl are added to Yeast Malt HiVeg™ Broth to eliminate moulds and thus permits enumeration of yeasts from mixed population.

Wickerham suggested the use of Yeast Malt Broth as an enrichment medium for yeasts by adding a layer of sterile paraffin oil (about 1 cm) on the surface of inoculated broth. After the growth occurs it should be streaked on YM Agar (MV424) to obtain isolated colonies of fermentative species. To isolate fermentative as well as oxidative strains, acidified YM Broth is placed on a rotary shaker for 1 or 2 days which favours development of yeast cells while the sporulation of moulds is prevented and yeasts can be isolated by streaking on Yeast Malt HiVeg™ Agar (YM HiVeg™ Agar) (MV424). Yeast Malt HiVeg™ Broth (YM HiVeg™ Broth) is prepared by using vegetable peptones in place of animal based peptones which make the media free of BSE/TSE risks.

HiVeg™ peptone serves as a source of carbon, nitrogen and essential nutrients. Yeast extract supplies vitamin B complex nutrients and other growth factors. Malt extract serves as an additional source of carbon. Dextrose is the carbohydrate and energy source. To increase the selectivity, the media can be acidified by the addition of sterile 10% HCl, tartaric acid or 10% citric acid. Alternatively antibiotics (penicillin 20U/ml or streptomycin to a final concentration of 40mcg/ml) can be added. Acidified medium should not be reheated.

Type of specimen

Food and dairy samples

Specimen Collection and Handling:

For food and dairy samples, follow appropriate techniques for sample collection and processing as per guidelines (4,5,7). After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions :

Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling clinical specimens. Safety guidelines may be referred in individual safety data sheets.

Limitations :

1. Further biochemical tests must be carried out for confirmation.
2. Acidified medium should not be reheated.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the expiry period when stored at recommended temperature.

Quality Control

Appearance

Cream to beige homogeneous free flowing powder

Colour and Clarity of prepared medium

Light amber coloured clear solution in tubes.

pH

6.00-6.40

Cultural Response

Cultural characteristics observed after an incubation at 25-30°C for 40-72 hours.

Organism	Growth at pH	
	3.4	6.2
<i>Aspergillus brasiliensis</i> ATCC 16404 (00053*)	good-luxuriant	good-luxuriant
<i>Candida albicans</i> ATCC 10231 (00054*)	good-luxuriant	good-luxuriant
<i>Escherichia coli</i> ATCC 25922 (00013*)	inhibited	good-luxuriant
<i>Lactobacillus leichmannii</i> ATCC 4797	poor	good-luxuriant
<i>Saccharomyces cerevisiae</i> ATCC 9763 (00058*)	good-luxuriant	good-luxuriant
<i>Lactobacillus casei</i> ATCC 9595	poor	good-luxuriant

* Corresponding WDCM numbers

Storage and Shelf Life

Store between 10-30°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition. Seal the container tightly after use. Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with sample must be decontaminated and disposed of in accordance with current laboratory techniques (6,7).

Reference

1. Wickerham, 1939, J. Tropical Med. Hyg., 42:176.
2. Wickerham, 1951, U.S. Dept. Agric. Tech. Bull. No. 1029.
3. American Public Health Association, Standard Methods for the Examination of Dairy Products, 1978, 14th Ed., Washington D.C.
4. Salfinger Y., and Tortorello M.L. Fifth (Ed.), 2015, Compendium of Methods for the Microbiological Examination of Foods, American Public Health Association, Washington, D.C.
5. Wehr H. M. and Frank J. H., 2004, Standard Methods for the Microbiological Examination of Dairy Products, 17th Ed., APHA Inc., Washington, D.C.
6. Isenberg, H.D. Clinical Microbiology Procedures Handbook 2nd Edition.
7. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.

Revision :04/2024

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory,diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.