

Spirit Blue HiVeg™ Agar

MV445

Intended Use:

Recommended for detection and enumeration of lipolytic microorganisms.

Composition**

Ingredients	g / L
HiVeg™ hydrolysate	10.000
Yeast extract	5.000
Spirit blue	0.150
Agar	17.000
Final pH (at 25°C)	6.8±0.2

**Formula adjusted, standardized to suit performance parameters

Directions

Suspend 32.15 grams in 1000 ml purified/distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C and add 30 ml lipase substrate slowly while agitating to obtain an even distribution. Mix well and pour into sterile Petri plates.

Principle And Interpretation

Lipids, including fats and oils, are highly reduced. When a lipid is catabolized, it has the potential to yield more pairs of electrons per gram, and thus more energy, than either carbohydrates or proteins (1). This process is brought about by the enzyme lipase, and the organisms possessing the enzyme lipase are called lipolytic organisms. Growth of lipase-producing microorganisms can contribute to flavour defects in milk and high fat dairy products. Some of the free fatty acids released by the action of lipolytic enzymes have a low flavour threshold and can impart a rancid flavour at low concentrations. Spirit Blue Agar is prepared according to the formulation of Starr (2) is recommended by APHA (3) for detection and enumeration of lipolytic microorganisms. It is a basal medium to which lipoidal substrate is added for the detection, enumeration and study of lipolytic microorganisms. Formulations in practice before Starr which included dyes as indicators of lipolysis were sometimes inhibitory to the microorganisms. Starr showed spirit blue to be inert and an ideal indicator of lipolysis, visualized as clear halos around colonies. Spirit Blue HiVeg™ Agar is prepared by using vegetable peptones in place of animal based peptones which make the media free of BSE/TSE risks. HiVeg™ hydrolysate and yeast extract in the medium are sources of carbon, nitrogen, vitamins and minerals. Spirit blue is a dye which acts as an indicator of lipolysis. The lipase reagents recommended as the lipid source are cotton seed meal, cream, olive oil etc. A satisfactory emulsion can be prepared by dissolving 10 gram acacia or 1 ml polysorbate 80 in 400 ml warm distilled water, adding 100 ml cotton seed or olive oil and agitating vigorously to emulsify.

Type of specimen

Isolated Microorganism from dairy samples

Specimen Collection and Handling:

Prepare 1:10 or other suitable dilution of the product to be tested. Spread 0.1 ml of the desired dilutions over the surface of the medium. Incubate at 35-37°C for 24-48 hours. Colonies of lipolytic organisms develop a clear zone and /or a deep blue colour around and under each colony (3,4).

Warning and Precautions :

Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheets.

Limitations :

1. For proper lipase activity, it is recommended to use glass plates instead of disposable plastic plates.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the expiry period when stored at recommended temperature.

Quality Control

Appearance

Cream to greenish yellow homogeneous free flowing powder

Gelling

Firm, comparable with 1.7% Agar gel.

Colour and Clarity of prepared medium

Basal medium yields blue coloured, clear to slightly opalescent gel. With addition of lipase substrate, lavender coloured slightly opalescent gel forms in Petri plates

Reaction

Reaction of 3.22% w/v aqueous solution at 25°C. pH : 6.8±0.2

pH

6.60-7.00

Cultural Response

Cultural characteristics observed after an incubation at 35-37°C for 48-72 hours with added Lipase substrate

Organism	Inoculum (CFU)	Growth	Lipase activity
<i>Proteus mirabilis</i> ATCC 25933	50-100	luxuriant	negative, absence of zone around colony
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> ATCC 25923 (00034*)	50-100	luxuriant	positive reaction, clear zone around colony
<i>Staphylococcus epidermidis</i> ATCC 12228 (00036*)	50-100	luxuriant	positive reaction, clear zone around colony

Key : *Corresponding WDCM numbers.

Storage and Shelf Life

Store between 10-30°C in a tightly closed container and the prepared medium at 20-30°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition. Seal the container tightly after use. Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with sample must be decontaminated and disposed of in accordance with current laboratory techniques (5,6).

Reference

1. Norton C. F., 1986, Microbiology, 2nd Ed., Addison-Wesley Publishing Company.
2. Starr, 1941, Science, 93:333.
3. Wehr H. M. and Frank J. H., 2004, Standard Methods for the Microbiological Examination of Dairy Products, 17th Ed., APHA Inc., Washington, D.C.
4. American Public Health Association, Standard Methods for the Examination of Dairy Products, 1978, 14th Ed., Washington D.C.
5. Isenberg, H.D. Clinical Microbiology Procedures Handbook 2nd Edition.
6. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.

Revision : 04/2024

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.