

Loeffler HiVeg™ Medium Base

MV537

Intended Use:

Recommended for cultivation of *Corynebacterium diphtheriae* from **clinical** specimens and in pure cultures, detection of chromogenesis, proteolysis and the production of ascospores.

Composition**

Ingredients	g / L
HiVeg™ special peptone	2.500
HiVeg™ extract	2.500
Sodium chloride	1.250
Dextrose (Glucose)	2.500
Final pH (at 25°C)	7.3±0.2

**Formula adjusted, standardized to suit performance parameters

Directions

Suspend 8.75 grams in 250 ml purified / distilled water. Dissolve the medium completely and sterilize by autoclaving at 10 lbs pressure (115°C) for 20 minutes. Cool to 50-55°C and aseptically add 750 ml of sterile horse serum (RM1239). Mix well and aseptically dispense into sterile tubes. Sterilize the medium by inspissation at 80-85°C for 2 hours in free flowing steam for atleast 3 consecutive days.

Principle And Interpretation

Loeffler Medium was originally devised by Loeffler (1) and was further modified by Perry and Petran (2) and Buck (3). Loeffler HiVeg™ Medium Base enhances primary and secondary isolation and cultivation of fastidious pathogenic microorganisms especially from nose and throat. It also restores virulence and other identifying properties (microscopic and colonial) after they have been lost due to prolonged incubation or repeated subculturing. Loeffler HiVeg™ Medium Base is prepared by completely replacing animal based peptone with vegetable peptones to avoid BSE/TSE risks associate with animal peptones. The high serum content helps in determining proteolytic activity of organisms. It is also used for demonstration of pigmentation and ascospores. HiVeg™ special peptone, HiVeg™ extract provide essential growth nutrients. Dextrose is the source of fermentable carbohydrate and energy. Rub the swabs directly over the surface of medium and after incubation; prepare the smears from surface of slope. For proteolysis testing, inoculate slant and prior to incubation, flood the slant with Brewer Thioglycollate HiVeg™ Medium (MV019). Incubation should be carried out for 3-4 days or much longer for appearance of proteolysis. Loeffler HiVeg™ Medium Base should be used in parallel with Serum Tellurite Agar for selective isolation of *Corynebacteria* (4).

Type of specimen

Clinical samples- a nasopharyngeal and throat swabs

Specimen Collection and Handling:

After use, contaminated materials must be sterilized by autoclaving before discarding.

Warning and Precautions:

Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheets.

Limitations :

1. For optimum recovery of *C. diphtheriae*, a nasopharyngeal and throat swab should be obtained while collection of specimen.
2. Further biochemical and serological tests need to be carried out for confirmation.

Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the expiry period when stored at recommended temperature.

Quality Control

Appearance

Cream to yellow homogeneous free flowing powder

Colour and Clarity of prepared medium

Basal: Light amber coloured clear to slightly opalescent solution. After addition of horse serum: Off-white coloured opaque solution

Reaction

Reaction of 3.52% w/v aqueous solution of base at 25°C. pH : 7.3±0.2

pH

7.10-7.50

Cultural Response

Cultural characteristics observed with added 750ml horse serum, after an incubation at 35-37°C for 24-48 hours.

Organism	Inoculum (CFU)	Growth
<i>Corynebacterium diphtheriae</i> ATCC 11913	50-100	fair-good
<i>Pseudomonas aeruginosa</i> ATCC 10145 (00024*)	50-100	good (green colonies with proteolysis)
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> ATCC 25923 (00034*)	50-100	good (yellow to gold colonies)

Key : *Corresponding WDCM numbers.

Storage and Shelf Life

Store between 10-30°C in a tightly closed container and the prepared medium at 2 - 8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition Seal the container tightly after use. Product performance is best if used within stated expiry period.

Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with sample must be decontaminated and disposed of in accordance with current laboratory techniques (5,6).

Reference

1. Loeffler F., 1887, Zentralb. Bakteriologie Parasitenkd., 2:102.
2. Perry and Petran, 1939, J. Lab. Clin. Med., 25:71.
3. Buck, 1949, J. Lab. Clin. Med., 34:582.
4. MacFaddin J. F., 1985, Media for Isolation-Cultivation-Maintenance of Medical Bacteria, Vol. I, Williams and Wilkins, Baltimore.
5. Isenberg, H.D. Clinical Microbiology Procedures Handbook 2nd Edition.
6. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.

Revision : 03/2024

Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.