

## MYP HiVeg™ Agar Base (Phenol Red Polymyxin HiVeg™ Agar Base)

MV636

### Intended Use:

Recommended for isolation and identification of *Bacillus* species and pathogenic Staphylococci.

### Composition\*\*

Ingredients	g / L
HiVeg™ peptone	10.000
HiVeg™ extract No. 1	1.000
D-Mannitol	10.000
Sodium chloride	10.000
Phenol red	0.025
Agar	15.000
Final pH ( at 25°C)	7.1±0.2

\*\*Formula adjusted, standardized to suit performance parameters

### Directions

Suspend 23.01 grams in 450 ml purified/distilled water. Heat to boiling to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45-50°C. Aseptically add rehydrated contents of one vial of sterile PolyB Selective Supplement (FD003) solution to a final concentration of 100 units per ml and 50 ml sterile Egg Yolk Emulsion (FD045). Mix well and pour into sterile Petri plates.

### Principle And Interpretation

*Bacillus cereus* is ubiquitously present in soil, vegetation water and dust. It has been isolated from a large variety of foods, including vegetables, meat, cereals, pasteurized fresh milk and powdered milk (1,2,3) and processed foods. Under favourable conditions, the organism multiplies and causes gastrointestinal illness (4). It is implicated in two different forms of food poisoning; an emetic illness and a diarrhoeal illness. The emetic illness is mediated by a highly stable toxin that survives high temperature, exposure to trypsin, pepsin and pH extremes. The diarrhoeal illness is mediated by a heat and acid labile enterotoxin.

Lecithinase activity is the key reaction in the differential identification of *B.cereus*, the most commonly encountered and important species in clinical laboratories, from the majority of the other *Bacillus* species. If unknown isolate produces lecithinase, *Bacillus cereus* can be presumptively identified by also observing colonial morphology, hemolytic reactivity and motility tests. Mossel et al (5) formulated Mannitol-Egg Yolk-Polymyxin (MYP) Agar, which is recommended by APHA to isolate and enumerate *B.cereus* from foods (3,4,6,7). When present in large numbers in certain foodstuffs, *B.cereus* can produce metabolites responsible for the clinical symptoms of food poisoning (5). This medium differentiates *B.cereus* from other bacteria based on the basis of lecithinase activity, mannitol fermentation and resistance to polymyxin (FD003) (3,8). MYP HiVeg™ Agar Base is prepared by completely replacing animal based peptone with vegetable peptones to avoid BSE/TSE risks associate with animal peptones. MYP HiVeg™ Agar Base contains HiVeg™ peptone and HiVeg™ extract No. 1, which provide nitrogen source. Mannitol fermentation can be detected by phenol red, which yields yellow colour to the mannitol fermenting colonies due to acid production. Added egg yolk emulsion helps in differentiation of lecithinase producing colonies, which are surrounded by a zone of white precipitate. Addition of Polymyxin B Sulphate (FD003) helps to restrict growth of gram-negative bacteria such as *Escherichia coli* and *Pseudomonas aeruginosa*. These differentiating media allow differentiation of *B.cereus* from other *Bacillus* species by its inability to ferment mannitol and poor sporulation. *B.cereus* dissimilates egg yolk and gives rise to typical bacilli form colonies.

### Type of specimen

Food samples; Water samples

### Specimen Collection and Handling:

For food samples, follow appropriate techniques for sample collection and processing as per guidelines (9). For water samples, follow appropriate techniques for sample collection, processing as per guidelines and local standards.(10) After use, contaminated materials must be sterilized by autoclaving before discarding.

## Warning and Precautions :

Read the label before opening the container. Wear protective gloves/protective clothing/eye protection/face protection. Follow good microbiological lab practices while handling specimens and culture. Standard precautions as per established guidelines should be followed while handling specimens. Safety guidelines may be referred in individual safety data sheets.

## Limitations :

1. If unknown isolate produces lecithinase, *Bacillus cereus* can be presumptively identified by also observing colonial morphology, hemolytic reactivity and motility tests.

## Performance and Evaluation

Performance of the medium is expected when used as per the direction on the label within the expiry period when stored at recommended temperature.

## Quality Control

### Appearance

Light yellow to light pink homogeneous free flowing powder

### Gelling

Firm, comparable with 1.5% Agar gel

### Colour and Clarity of prepared medium

Basal medium :Red coloured clear to slightly opalescent gel. After Addition of Egg Yolk Emulsion (FD045) : Light orange coloured opaque gel forms in Petri plates

### Reaction

Reaction of 4.6% w/v aqueous solution at 25°C. pH : 7.1±0.2

### pH

6.90-7.30

### Cultural Response

Cultural characteristics observed with added Egg Yolk Emulsion (FD045)and PolyB Selective Supplement(FD003) after an incubation at 32°C for 18-40 hours.

Organism	Inoculum (CFU)	Growth	Recovery	Colour of colony	Lecithinase activity
<i>Bacillus cereus</i> ATCC 10876	50-100	luxuriant	≥50%	red	positive, opaque zone around the colony
** <i>Bacillus spizizenii</i> ATCC 6633 (00003*)	50-100	luxuriant	≥50%	yellow	negative
<i>Escherichia coli</i> ATCC 25922 (00013*)	50-100	none-poor	≤10%		
<i>Proteus mirabilis</i> ATCC 25933	50-100	luxuriant	≥50%	red	negative
<i>Pseudomonas aeruginosa</i> ATCC 27853 (00025*)	50-100	none-poor	≤10%		
<i>Staphylococcus aureus</i> subsp. <i>aureus</i> ATCC 25923 (00034*)	50-100	luxuriant	≥50%	yellow	positive, opaque zone around the colony

Key : \*Corresponding WDCM numbers and \*\*Formerly known as *Bacillus subtilis* subsp. *spizizenii*.

## Storage and Shelf Life

Store between 10-30°C in a tightly closed container and the prepared medium at 2-8°C. Use before expiry date on the label. On opening, product should be properly stored dry, after tightly capping the bottle in order to prevent lump formation due to the hygroscopic nature of the product. Improper storage of the product may lead to lump formation. Store in dry ventilated area protected from extremes of temperature and sources of ignition. Seal the container tightly after use. Product performance is best if used within stated expiry period.

## Disposal

User must ensure safe disposal by autoclaving and/or incineration of used or unusable preparations of this product. Follow established laboratory procedures in disposing of infectious materials and material that comes into contact with sample must be decontaminated and disposed of in accordance with current laboratory techniques (11,12).

## Reference

1. Bergdoll M. S., 1981, Clin. Microbiol. Newsletter 3: 85-87.
2. Centers for Disease Control: Bacillus cereus- Maine, MMWR, 35: 408-410, 1986.
3. Downes F. P. and Ito K., (Eds.), 2001, Compendium of Methods for the Microbiological Examination of Foods, 4th Ed., APHA, Washington, D.C.
4. Donovan K. O., 1958, J. Appl. Bacteriol., 21:100.
5. Mossel D. A. A., Koopman M. J. and Jongerium E., 1967, Appl. Microbiol, 15:650.
6. Colimer A. R., 1948, J. Bacteriol., 55:777.
7. Rhodehamel J. and Harmon S. M., 1995, FDA Bacteriological Analytical Manual, 8th Ed., AOAC International, Gaithersburg, Md.
8. Nygren B., 1962, Acta Path. Microbiol. Scand., 56: Suppl. 1.
9. Salfinger Y., and Tortorello M.L., 2015, Compendium of Methods for the Microbiological Examination of Foods, 5th Ed., American Public Health Association, Washington, D.C.
10. Baird R.B., Eaton A.D., and Rice E.W., (Eds.), 2015, Standard Methods for the Examination of Water and Wastewater, 23rd ed., APHA, Washington, D.C.
11. Isenberg, H.D. Clinical Microbiology Procedures Handbook 2nd Edition.
12. Jorgensen, J.H., Pfaller, M.A., Carroll, K.C., Funke, G., Landry, M.L., Richter, S.S and Warnock., D.W. (2015) Manual of Clinical Microbiology, 11th Edition. Vol. 1.

Revision : 03/2024

### Disclaimer :

User must ensure suitability of the product(s) in their application prior to use. Products conform solely to the information contained in this and other related HiMedia™ publications. The information contained in this publication is based on our research and development work and is to the best of our knowledge true and accurate. HiMedia™ Laboratories Pvt Ltd reserves the right to make changes to specifications and information related to the products at any time. Products are not intended for human or animal or therapeutic use but for laboratory, diagnostic, research or further manufacturing use only, unless otherwise specified. Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.