

---

## Potassium phosphate dibasic anhydrous, Plant Culture Tested PCT0018

---

### Product Identifier

---

CAS No.	:	7758-11-4
EC No.	:	231-834-5
Molecular Formula	:	K <sub>2</sub> HPO <sub>4</sub>
Molecular Weight	:	174.18
HS Code	:	2835 24 00
Storage	:	Below 30°C
Shelf life	:	4 years

### Technical Specification

---

Appearance	:	Colourless or white to off-white hygroscopic crystals or powder or granules or solid
Solubility	:	1000 mg soluble in 1 mL of water
pH (5% in water at 25°C)	:	8.50 - 9.60
Loss on drying (at 105°C, constant weight)	:	<= 1.0%
Assay (NaOH/HCl Titration, on dry basis)	:	99.00 - 101.00%
Plant Tissue Culture Test	:	Passes

### Safety Information

---

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1