

---

## EDTA ferric monosodium salt, Plant Culture Tested

**PCT0106**

### Product Identifier

---

CAS No.	:	15708-41-5
EC No.	:	239-802-2
Molecular Formula	:	C <sub>10</sub> H <sub>12</sub> N <sub>2</sub> O <sub>8</sub> FeNa
Molecular Weight	:	367.05
HS Code	:	2922 49 90
Storage	:	Below 30°C
Shelf life	:	4 years

### Technical Specification

---

Appearance	:	Light yellow to yellowish green to brown powder
Solubility	:	Clear solution at 33.3 mg in 1 mL of water
Assay (Fe, Iodometry)	:	12.00 - 14.00%
Plant Tissue Culture Test	:	Passes

### Safety Information

---

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	AH4900000
WGK	:	2