

## Ferric ammonium citrate, Brown, Plant Culture Tested

PCT0107

### Product Identifier

CAS No.	:	1185-57-5
EC No.	:	214-686-6
Synonym	:	Ammonium ferric citrate, Brown
HS Code	:	2918 15 50
Storage	:	Below 30°C
Shelf life	:	4 years

### Technical Specification

Appearance	:	Yellow to reddish-brown crystals or powder or granules
Solubility	:	100 mg soluble in 1 mL of water
Ferric citrate	:	Passes test
Assay (% Fe, Iodometry)	:	16.50 - 22.50%
Plant Tissue Culture Test	:	Passes

### Safety Information

Hazard Pictogram(s)



Signal Word	:	Warning
Hazard Statement(s)	:	H315- H319- H335
Precautionary Statement(s)	:	P261- P280- P312- P321- P403 + P233- P501
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	GE7540000
WGK	:	1