

Molybdenum trioxide Plant Culture Tested

PCT0114


Product Identifier

CAS No.	:	1313-27-5
EC No.	:	215-204-7
Molecular Formula	:	MoO ₃
Molecular Weight	:	143.94
HS Code	:	2825 70 10
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to slightly yellow to green or grey crystals or powder or solid
Solubility	:	33.3 mg soluble in 1 mL of dilute sodium hydroxide solution
Assay (Lead nitrate Titration)	:	min. 99.00%
Plant Tissue Culture Test	:	Passes

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H319- H335- H351
Precautionary Statement(s)	:	P261- P281- P305+P351+P338
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	QA4725000
WGK	:	3