
Ammonium molybdate tetrahydrate, Plant Culture Tested PCT0119


Product Identifier

CAS No.	:	12054-85-2
EC No.	:	234-320-9
Molecular Formula	:	$(\text{NH}_4)_6\text{Mo}_7\text{O}_{24}\cdot 4\text{H}_2\text{O}$
Molecular Weight	:	1235.86
HS Code	:	2841 70 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White crystals or powder, sometimes with a slight green or yellow tint
Solubility	:	33.3 mg soluble in 1 mL of water
Assay (Lead nitrate Titration)	:	99.00 - 102.00%
Plant Tissue Culture Test	:	Passes

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H302- H315- H319- H335
Precautionary Statement(s)	:	P261- P280- P312- P321- P403 + P233- P501
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	QA5076000
WGK	:	1