
Sodium thiosulphate pentahydrate, Plant Culture Tested

PCT0122

Product Identifier

CAS No.	:	10102-17-7
EC No.	:	231-867-5
Molecular Formula	:	Na ₂ S ₂ O ₃ .5H ₂ O
Molecular Weight	:	248.18
HS Code	:	2832 30 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Colourless or white crystals or powder or pellets or solid
Solubility	:	1000 mg soluble in 1 mL of water
pH (5% in water at 25°C)	:	5.50 - 8.40
Assay (Iodimetry)	:	min. 99.50%
Plant Tissue Culture Test	:	Passes

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1