

---

## Choline chloride, Plant Culture Tested

**PCT0203**

---

### Product Identifier

CAS No.	:	67-48-1
EC No.	:	200-655-4
Molecular Formula	:	C <sub>5</sub> H <sub>14</sub> ClNO
Molecular Weight	:	139.62
HS Code	:	2923 10 00
Storage	:	Below 30°C
Shelf life	:	3 years

---

### Technical Specification

Appearance	:	White deliquescent crystals or powder
Solubility	:	Clear solution at 33.3 mg in 1 mL of water
FTIR	:	Matches with the standard pattern
Water (K.F.)	:	<= 5.00%
Assay (AT/NT, anhydrous basis)	:	99.00 - 102.00%
Plant Tissue Culture Tested	:	Passes

---

### Safety Information

Hazard Pictogram(s) :



Signal Word	:	Warning
Hazard Statement(s)	:	H315- H319- H335
Precautionary Statement(s)	:	P261- P280- P312- P403 + P233- P501
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	KH2975000
WGK	:	1