

PCT0401 Apple Powder

Product Number	Packing
PCT0401	: 250G

Product Information

Product Code	: PCT0401
Product Name	: Apple Powder
Shelf Life	: 4 years

Technical Specification

Appearance	: Cream to dark cream coloured powder
Reaction	: Reaction of 5.0% w/v aqueous solution at 25°C
pH	: min. 4.00 max. 5.00
Cultural response	: Cultures conditions - Incubation period (5wks), Relative humidity (60±2%), Temperature (25±2°C), Photoperiod Day: Night in hours (16:8)
Shoot culture	: No structural deformity observed, actively growing shoots, no toxicity to shoots
Moisture content	: 6.0% (max)
Total ash	: 6.0% (max)
Acidity (as Malic acid)	: 3.0% (max)

Risk and Safety Information

Storage Temperature(°C)	: Store below 30°C
-------------------------	--------------------

Transport Information

ADR/RID	: Not Dangerous Goods
IMDG	: Not Dangerous Goods
IATA	: Not Dangerous Goods