
Banana powder, Plant Culture Tested

PCT0402

Product Identifier

HS Code	:	1106 30 90
Storage	:	Below 30°C
Shelf life	:	2 years

Technical Specification

Appearance	:	White to tan homogeneous free flowing powder
Solubility	:	Insoluble in water
pH (3% in water at 25°C)	:	4.00 - 6.00
Water (K.F.)	:	<= 8.00%
Sulphated ash	:	<= 3.00%
Plant Tissue Culture Test	:	Passes

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	3