



## Product Information

### PCT0409 Potato Powder

---

Product Number	Packing
----------------	---------

---

PCT0409	: 100G
---------	--------

---

#### Product Information

---

Product Code	: PCT0409
Product Name	: Potato Powder
Shelf Life	: 4 years

---

#### Technical Specification

---

Appearance	: Light brown coloured, free flowing coarse powder
Solubility	: Insoluble in cold water & soluble in water at 45°C
Reaction	: Reaction of 1.0% w/v aqueous solution at 25°C
pH	: min. 4.90 max. 5.90
Cultural response	: Cultures conditions - Incubation period (5wks), Relative humidity (60±2%), Temperature (25±2°C), Photoperiod Day: Night in hours (16:8)
Shoot culture	: No structural deformity observed, actively growing shoots, no toxicity to shoots
Moisture % (KF)	: <= 8.0%
Ash %	: <= 3.50%
Sulphate (SO <sub>4</sub> )	: <= 0.003%
alpha amino -N	: <= 1.70%

---

#### Risk and Safety Information

---

Storage Temperature(°C)	: store below 30°C
-------------------------	--------------------

---

#### Transport Information

---

ADR/RID	: Not Dangerous Goods
IMDG	: Not Dangerous Goods
IATA	: Not Dangerous Goods