

Citric acid anhydrous, Plant Culture Tested

PCT0501


Product Identifier

CAS No.	:	77-92-9
EC No.	:	201-069-1
Molecular Formula	:	C ₆ H ₈ O ₇
Molecular Weight	:	192.12
HS Code	:	2918 14 00
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Colourless or white to almost white hygroscopic crystals or powder or solid or granules
Solubility	:	100 mg soluble 1 mL of water
Water (K.F.)	:	<= 1.00%
Assay (NaOH Titration, on anhydrous basis)	:	>= 99.70%
Plant Tissue Culture Test	:	Passes

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H319- H335
Precautionary Statement(s)	:	P261- P271- P280- P304 + P340- P305+P351+P338- P403 + P233- P501
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	GE7350000
WGK	:	1