
Maltodextrin, Plant Culture Tested

PCT0611

Product Identifier

CAS No.	:	9050-36-6
EC No.	:	232-940-4
HS Code	:	3505 10 10
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White to off white, slightly hygroscopic powder or granules
Solubility	:	100 mg soluble in 1 mL of water giving clear to hazy solution
Identification	:	Complies (To 5 mL of hot alkaline cupric tartrate TS add a few drops of a 1 in 10 solution of the sample. A red precipitate of cuprous oxide form)
Reducing sugar (Dextrose equivalent)	:	<= 20.00%
Plant Tissue Culture Test	:	Passes

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	3