
trans-Zeatin riboside, Plant Culture Tested

PCT0814

Product Identifier

CAS No.	:	6025-53-2
Molecular Formula	:	C ₁₅ H ₂₁ N ₅ O ₅
Molecular Weight	:	351.36
HS Code	:	2934 99 00
Storage	:	2 - 8°C
Shelf life	:	2 years

Technical Specification

Appearance	:	White to off-white to tan crystals or powder
Solubility	:	33.3 mg soluble in 1 mL of sodium hydroxide
Melting range	:	180 - 182°C
Assay (HPLC)	:	min. 98.00%
Plant Tissue Culture Test	:	Passes

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	3