
meta-Topolin, Plant Culture Tested

PCT0827

Product Identifier

CAS No.	:	75737-38-1
Molecular Formula	:	C ₁₂ H ₁₀ N ₅ OH
Molecular Weight	:	241.25
HS Code	:	2933 59 90
Storage	:	2 - 8°C
Shelf life	:	3 years

Technical Specification

Appearance	:	White to yellow crystals or powder
Solubility	:	33.3 mg soluble in 1 mL of 1N sodium hydroxide
Plant Tissue Culture Test	:	Passes

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available