
Ancymidol, Plant Culture Tested

PCT0837

Product Identifier

CAS No.	:	12771-68-5
EC No.	:	235-814-7
Molecular Formula	:	C ₁₅ H ₁₆ N ₂ O ₂
Molecular Weight	:	256.30
HS Code	:	2933 59 90
Storage	:	2 - 8°C
Shelf life	:	3 years

Technical Specification

Appearance	:	White to light yellow coarse crystals or powder or solid
Solubility	:	33.3 mg soluble in 1 mL of DMSO
Assay (HPLC/GC)	:	min. 98.00%
Plant Tissue Culture Test	:	Passes

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	UV9280000
WGK	:	1