

---

## Flurprimidol, Plant Culture Tested

PCT0844

### Product Identifier

---

CAS No.	:	56425-91-3
Molecular Formula	:	C <sub>15</sub> H <sub>15</sub> F <sub>3</sub> N <sub>2</sub> O <sub>2</sub>
Molecular Weight	:	312.29
HS Code	:	2933 59 90
Storage	:	Below 30°C
Shelf life	:	4 years


### Technical Specification

---

Appearance	:	Colourless or white to dark beige crystals or powder
Solubility	:	Clear solution at 33.3 mg in 1 mL of methanol
Assay (HPLC)	:	min. 98.00%
Plant Tissue Culture Test	:	Passes

### Safety Information

---

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H302- H315- H319
Precautionary Statement(s)	:	P270- P280- P301 + P312- P302 + P352- P305+P351+P338- P501
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1