

Product Information

PCT0902 CleriGarTM

Product Number		Packing
PCT0902	:	1KG
PCT0902	:	25KG
PCT0902	:	500G
PCT0902	:	5KG

Product Information

Product Code : PCT0902
Product Name : CleriGarTM
Shelf Life : 4 years

Technical Specification

Appearance : White to off-white coloured, homogeneous, free flowing powder

Solubility : 1.5% (w/v) soluble after boiling in distilled water

Clarity : A firm solid, clear to slight opalescent gel is formed at a concentration of

1.5% at 35 - 40°C

Reaction : Reaction of 1.5% w/v aqueous solution at 25 °C

pH : min. 5.50 max. 7.50

Cultural response : Cultures conditions - Incubation period (5wks), Relative humidity (60±2%),

Temperature (25±2°C), Photoperiod Day: Night in hours (16:8)

Shoot culture : No structural deformity observed, actively growing shoots, no toxicity to

shoots

Callus culture : No necrotic tissues, actively growing callus, no toxicity to callus

Gelling temperature : 35 - 40°C

Moisture : Not more than 20%

Risk and Safety Information

Storage Temperature(°C) : Store below 30°C

Transport Information

ADR/RID : Not Dangerous Goods IMDG : Not Dangerous Goods IATA : Not Dangerous Goods

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919 Fax:+91-022-6147 1920 Gram: STERILITY.email:info@himedialabs.com. Visit us at our website: www.himedialabs.com