



## Product Information

### PCT0903 CleriGel™ Plant Culture Tested

#### Product Number Packing

PCT0903	: 100G
PCT0903	: 1KG
PCT0903	: 5KG
PCT0903	: 25KG

#### Product Information

Product Code	: PCT0903
Product Name	: CleriGel™; Plant Culture Tested
Synonym	: Gellan gum
CAS No.	: 71010-52-1
EC No.	: 275-117-5
HS Code	: 3913 90 00
Shelf Life	: 4 years

#### Technical Specification

Appearance	: White to off-white coloured, homogeneous, free flowing powder
Solubility	: Freely soluble in water at temperatures above 85°C
Colour	: Colourless to off white
Clarity	: Clear to slight hazy
Reaction	: Reaction of 0.2% w/v aqueous solution at 25°C
pH	: min. 5.50 max. 7.50
Cultural response	: Cultures conditions - Incubation period (5wks), Relative humidity (60±2%), Temperature (25±2°C), Photoperiod Day: Night in hours (16:8)
Shoot culture	: No structural deformity observed, actively growing shoots, no toxicity to shoots
Callus culture	: No necrotic tissues, actively growing callus, no toxicity to callus
Loss on drying	: ≤ 15.0%

#### Risk and Safety Information

WGK	: 2
Storage Temperature(°C)	: Store below 30°C

#### Transport Information

Marine Pollutant	: No
ADR/RID	: Not Dangerous Goods
IMDG	: Not Dangerous Goods
IATA	: Not Dangerous Goods

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919  
Fax:+91-022-6147 1920 Gram : STERILITY . email : [info@himedialabs.com](mailto:info@himedialabs.com). Visit us at our website: [www.himedialabs.com](http://www.himedialabs.com)

The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such Information or with respect to any product, method or apparatus referred to herein