
Antimicrobial supplement Plant Culture Tested

PCT1106

Product Identifier

HS Code : 2941 90 90
Storage : On receipt store between 10-30°C
Shelf life : 3 years

Product Information

Product Name : Antimicrobial Supplement
Plant Culture Tested

Antimicrobial supplement is a biodegradable antimicrobial product. It has a wide range of activity and can be used against various microbes. It has been tested at 0.05-0.5% concentration on various plant species with various microbes.

Technical Specification

Appearance : Colourless to yellow liquid
Solubility : Miscible in water
pH : min 2.00 - max 4.00
Cultural response : Antibacterial and Antifungal activity is determined by using

Soyabean Caesin Agar incubated at 30-35°C for 18-24 hours for bacteria and Sabouraud Dextrose Agar incubated at 20-25°C for 3-5 days for fungi and yeasts.

Safety Information

Hazard Pictogram(s) :



Hazard Statement(s) : H315- H319- H335
Precautionary Statement(s) : P264- P280- P305+P351+P338- P332/P337 + P313
UN No. :
Class : -
Packing Group : -