
Nystatin, Plant Culture Tested

PCT1112

Product Identifier

CAS No.	:	1400-61-9
EC No.	:	215-749-0
Molecular Formula	:	C ₄₇ H ₇₅ NO ₁₇
Molecular Weight	:	926.09
HS Code	:	2941 90 90
Storage	:	2 - 8°C
Shelf life	:	3 years

Technical Specification

Appearance	:	Off-white to slightly brownish hygroscopic powder
Solubility	:	Clear solution at 100 mg in 1 mL of dimethylformamide
Potency	:	>= 4400 USP units/mg

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	RF5950000
WGK	:	1