

Timentin, Plant Culture Tested

PCT1113

Ticarcillin disodium salt : Potassium clavulanate mixture (15:1)


Product Identifier

CAS No.	:	86482-18-0
HS Code	:	2941 90 90
Storage	:	2 - 8°C
Shelf life	:	3 years

Technical Specification

Appearance	:	White to yellow powder.
Solubility	:	Clear solution at 33.3 mg in 1ml of water
pH (0.005% in water at 25°C)	:	5.50 - 7.50
Plant Tissue Culture Test	:	Passes

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H315- H319- H334- H335
Precautionary Statement(s)	:	P261- P280- P284- P304 + P340- P312- P501
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1