

Vancomycin hydrochloride, Plant Culture Tested

PCT1114


Product Identifier

CAS No.	:	1404-93-9
Molecular Formula	:	C ₆₆ H ₇₅ Cl ₂ N ₉ O ₂₄ .HCl
Molecular Weight	:	1485.71
HS Code	:	2941 90 90
Storage	:	2 - 8°C
Shelf life	:	3 years

Technical Specification

Appearance	:	White to brown crystals or powder
Solubility	:	100 mg soluble in 1 mL of water
FTIR	:	Matches with the standard pattern
Potency	:	>= 950 µg/mg Vancomycin
Assay (HPLC)	:	min. 80.00%
Plant Tissue Culture Test	:	Passes

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H317- H334
Precautionary Statement(s)	:	P261- P280- P284- P304 + P340- P362 + P364- P501
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	YW4380000
WGK	:	1