

---

## Amoxicillin:clavulanic acid potassium salt (5:1), Plant Culture Tested

**PCT1115**

---

### Product Identifier

Synonym	:	Augmentin
HS Code	:	3004 90 99
Storage	:	2 - 8°C
Shelf life	:	3 years

---

### Technical Specification

Appearance	:	White to yellow powder
Solubility	:	33.3 mg soluble in 1 mL of water
pH	:	8.00 - 10.00
Plant Tissue Culture Test	:	Passes

---

### Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1