

---

## Erythromycin, Plant Culture Tested

**PCT1123**

### Product Identifier

---

CAS No.	:	114-07-8
EC No.	:	204-040-1
Molecular Formula	:	C <sub>37</sub> H <sub>67</sub> NO <sub>13</sub>
Molecular Weight	:	733.93
HS Code	:	2941 50 00
Storage	:	Below 30°C
Shelf life	:	3 years

### Technical Specification

---

Appearance	:	Colourless or white to yellow hygroscopic crystals or powder
Solubility	:	Clear solution at 100 mg in 1 mL of alcohol
Potency	:	>= 850 µg/mg
Plant Tissue Culture Test	:	Passes

### Safety Information

---

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	KF4375000
WGK	:	1