
Colistin sulphate, Plant Culture Tested

PCT1142

Product Identifier

CAS No.	:	1264-72-8
EC No.	:	215-034-3
HS Code	:	2941 90 90
Storage	:	2 - 8°C
Shelf life	:	3 years

Technical Specification

Appearance	:	White to pale yellow crystals or powder
Solubility	:	Clear solution at 33.3 mg in 1 mL of water
pH (1% in water at 25°C)	:	4.00 - 7.00
Potency	:	>= 11500 units/mg
Plant Tissue Culture Test	:	Passes

Safety Information

Hazard Pictogram(s)



Signal Word	:	Danger
Hazard Statement(s)	:	H301
Precautionary Statement(s)	:	P264.- P270- P301 + 310- P330- P405- P501
UN No.	:	2811
Class	:	6.1
Packing Group	:	III
RTECS	:	TR1500000
WGK	:	1