
BCIP®, Plant Culture Tested

PCT1201

Product Identifier

CAS No.	:	102185-33-1
Molecular Formula	:	C ₈ H ₄ BrClNO ₄ PNa ₂
Molecular Weight	:	370.43
Synonym	:	5-Bromo-4 chloro-3-indolyl phosphate disodium salt
HS Code	:	2933 99 90
Storage	:	2 - 8°C
Shelf life	:	2 years

Technical Specification

Appearance	:	White to light beige or light blue crystals or powder
Solubility	:	33.3 mg soluble in 1 mL of water
Water (K.F.)	:	6.0 - 10.0%
Assay (NT/HPLC)	:	98.00 - 102.00%
Plant Tissue Culture Test	:	Passes

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1