

---

## X-Glucuronide CHA salt, Plant Culture Tested

**PCT1205**

### Product Identifier

---

CAS No.	:	114162-64-0
Molecular Formula	:	C <sub>20</sub> H <sub>26</sub> BrClN <sub>2</sub> O <sub>7</sub>
Molecular Weight	:	521.79
Synonym	:	5-Bromo-4-chloro-3-indolyl-β-D-glucuronide cyclohexylammonium salt
HS Code	:	2934 99 90
Storage	:	2 - 8°C
Shelf life	:	2 years

### Technical Specification

---

Appearance	:	White to off-white crystals or powder or solid
Solubility	:	33.3 mg soluble in 1 mL of dimethyl formamide
Specific rotation	:	-95.0° to -85.0° (c = 1% in Dimethylformamide at 25°C, anhydrous basis)
Water (K.F.)	:	<= 6.0%
Assay (HPLC/NT, anhydrous basis)	:	98.00 - 102.00%
Plant Tissue Culture Test	:	Passes

### Safety Information

---

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1