
Acetocarmine, Plant Culture Tested

PCT1304

Product Identifier

HS Code : 3203 00 90
Storage : Below 30°C
Shelf life : 2 years

Technical Specification

Appearance : Deep red to orangish red liquid.
Solubility : 33.3 mg soluble in 1 mL of water
Plant Tissue Culture Test : Passes

Safety Information

Hazard Pictogram(s) :



Signal Word : Danger
Hazard Statement(s) : H226- H314
Precautionary Statement(s) : P210- P280- P303+P361+P353- P305+P351+P338- P370
+ P378- P403 + P235
UN No. : 2790
Class : 8
Packing Group : III
RTECS : Not available
WGK : 1