



Product Information

PCT1308 Silver thiosulphate Plant Culture Tested

Product Number	Packing
PCT1308	: 1X100ML
PCT1308	: 5X100ML

Product Information

Product Code	: PCT1308
Product Name	: Silver thiosulphate; Plant Culture Tested
Shelf Life	: 6 months
Storage	: Store at 2-8 ^o C

Technical Specification

Colour	: Colourless
Clarity	: Clear solution
Culture response	: Cultures conditions - Incubation period (5wks), Relative humidity (60±2%), Temperature (25±2°C), Photoperiod Day: Night in hours (16:8)
Shoot culture	: No structural deformity observed, actively growing shoots, no toxicity to shoots
Callus culture	: No necrotic tissues, actively growing callus, no toxicity to callus.
Sterility	: No bacterial or fungal growth is observed after 14 days of incubation, as per USP specification.

A-516, Swastik Disha Business Park, Via Vadhani Indl. Est., LBS Marg, Mumbai-400 086, India. Phone:+91-022-6147 1919
Fax:+91-022-6147 1920 Gram : STERILITY . email : info@himedialabs.com. Visit us at our website: www.himedialabs.com

.....
The information contained herein is believed to be accurate and complete. However no warranty or guarantee whatsoever is made or is to be implied with respect to such information or with respect to any product, method or apparatus referred to herein