

Benzalkonium chloride, 50% (w/v) aqueous solution, Plant Culture Tested

PCT1322


Product Identifier

CAS No.	:	63449-41-2
HS Code	:	3402 41 00
Storage	:	Below 30°C
Shelf life	:	2 years

Technical Specification

Appearance	:	Colourless to pale yellow liquid
Solubility	:	1 mL miscible in 1 mL of water
Assay (AT)	:	48.00 - 55.00%
Plant Tissue Culture Test	:	Passes

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H302+H312- H314- H400
Precautionary Statement(s)	:	P260- P273- P280- P301+P330+P331- P303+P361+P353- P305+P351+P338- P310- P501
UN No.	:	3265
Class	:	8
Packing Group	:	II
RTECS	:	BO3151000
WGK	:	2