

Sodium hydroxide pellets, Plant Culture Tested

PCT1325

Product Identifier

CAS No.	:	1310-73-2
EC No.	:	215-185-5
Molecular Formula	:	NaOH
Molecular Weight	:	40.00
HS Code	:	2815 11 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White hygroscopic pellets
Solubility	:	100 mg soluble in 1 mL of water
Carbonate (Na ₂ CO ₃)	:	<= 2.00%
Assay (HCl Titration)	:	97.00 - 100.50%
Plant Tissue Culture Test	:	Passes

Safety Information

Hazard Pictogram(s)



Signal Word	:	Danger
Hazard Statement(s)	:	H314
Precautionary Statement(s)	:	P260- P280- P305+P351+P338- P310
UN No.	:	1823
Class	:	8
Packing Group	:	II
RTECS	:	WB4900000
WGK	:	1