

Mercuric chloride, Plant Culture Tested

PCT1512


Product Identifier

CAS No.	:	7487-94-7
EC No.	:	231-299-8
EC Index No	:	080-010-00-X
Molecular Formula	:	HgCl ₂
Molecular Weight	:	271.5
HS Code	:	2852 90 00
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Colourless to white crystals or powder or solid
Solubility	:	33.3 mg soluble in 1 mL of water
Assay (EDTA Titration)	:	99.50 - 102.00%
Plant Tissue Culture Test	:	Passes

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Danger
Hazard Statement(s)	:	H300- H314- H341- H361f- H372- H410
Precautionary Statement(s)	:	P260- P273- P280- P301 + 310- P303+P361+P353- P305+P351+P338
UN No.	:	1624
Class	:	6.1
Packing Group	:	II
RTECS	:	OV9100000
WGK	:	3