

## Methylene blue trihydrate, Plant Culture Tested

**PCT1517**


### Product Identifier

CAS No.	:	7220-79-3
EC No.	:	200-515-2
C.I. No.	:	52015
Molecular Formula	:	C <sub>16</sub> H <sub>18</sub> Cl N <sub>3</sub> S.3H <sub>2</sub> O
Molecular Weight	:	373.90
HS Code	:	3204 13 51
Storage	:	Below 30°C
Shelf life	:	4 years

### Technical Specification

Appearance	:	Green to dark green crystals or powder with metallic lustre
Solubility	:	Clear blue solution at 20 ppm in water
Absorption maxima	:	658 - 668 nm (20 ppm in water)
Loss on drying (at 105°C, 2 hr)	:	8.00 - 18.00%
Assay (AT/UV, on dry basis)	:	min. 99.50%
Plant Tissue Culture Test	:	Passes

### Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H302
Precautionary Statement(s)	:	P264- P270- P301 + P312- P330- P501
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	SP5740000
WGK	:	1