
Pectinase, From *Aspergillus niger*, Plant Culture Tested

PCT1519

Product Identifier

CAS No.	:	9032-75-1
EC No.	:	232-885-6
HS Code	:	2922 50 90
Storage	:	2 - 8°C
Shelf life	:	3 years

Technical Specification

Appearance	:	White to slightly beige to pale amber solid or powder
Solubility	:	Clear solution at 33.3 mg in 1 mL of water
Loss on drying (at 105°C, 1 hr)	:	<= 10.00%
Activity	:	8000 - 12000 U/g
Plant Tissue Culture Test	:	Passes
Unit definition	:	One unit will liberate 1.0 micromole of galacturonic acid from polygalact-uronic acid per minute at pH 4.0 at 25°C

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available