
Sodium hydrogen carbonate, Plant Culture Tested

PCT1535

Product Identifier

CAS No.	:	144-55-8
EC No.	:	205-633-8
Molecular Formula	:	NaHCO ₃
Molecular Weight	:	84.01
Synonym	:	Sodium bicarbonate
HS Code	:	2836 30 00
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	White crystals or powder or solid.
Solubility	:	33.3 mg soluble in 1 mL of water
Assay (HCl Titration)	:	99.00 - 102.00%
Plant Tissue Culture Test	:	Passes

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	VZ0950000
WGK	:	1