

DAPI dihydrochloride Plant Culture Tested

PCT1539

Product Identifier

CAS No.	:	28718-90-3
EC No.	:	249-186-7
Molecular Formula	:	C ₁₆ H ₁₅ N ₅ .2HCl
Molecular Weight	:	350.25
Synonym	:	4',6-Diamidino-2-phenylindole dihydrochloride
HS Code	:	2933 99 00
Storage	:	2 - 8°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Yellow crystals or powder or solid
Solubility	:	Clear solution at 10 mg in 1 mL of water
Assay (HPLC)	:	min. 98.00%
Plant Tissue Culture Test	:	Passes

Safety Information

UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	Not available
WGK	:	1