

Dimethylsiloxane, ethylene oxide block copolymer Plant Culture Tested

PCT1554


Product Identifier

CAS No.	:	27306-78-1
Molecular Formula	:	C ₁₃ H ₃₄ O ₄ Si ₃
Molecular Weight	:	338.66
Synonym	:	Silwet® L-77
HS Code	:	2931 90 90
Storage	:	Below 30°C
Shelf life	:	4 years

Technical Specification

Appearance	:	Colourless to slight amber liquid.
Solubility	:	Dispersable in water
Specific Gravity	:	1.02 - 1.04
Plant Tissue Culture Test	:	Passes

Safety Information

Hazard Pictogram(s)	:	
Signal Word	:	Warning
Hazard Statement(s)	:	H302- H315- H319
Precautionary Statement(s)	:	P280- P305+P351+P338
UN No.	:	Not dangerous goods
Class	:	-
Packing Group	:	-
RTECS	:	OO9650000