

PGE0008 HiGroGA Powder

Product Number Packing

PGE0008 : 10G
PGE0008 : 25G

Product Information

Product Code : PGE0008
Product Name : HiGroGA Powder
Shelf Life : 2 years

Technical Specification

Appearance : White to pale yellow crystals or powder
Solubility : Soluble in 1N NaOH
Melting range : ~227°C
Assay (NaOH Titration/HPLC) : min. 90.00 %

Risk and Safety Information

R-Phrase(s) : 36
S-Phrase(s) : 26-36
Storage Temperature(°C) : Store below 30°C
Hazard Symbol(s)



Transport Information

Marine Pollutant : No
ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods