

PGE0009 HiGroBA Powder

Product Number Packing

PGE0009 : 5G
PGE0009 : 25G

Product Information

Product Code : PGE0009
Product Name : HiGroBA Powder
Shelf Life : 2 years

Technical Specification

Appearance : White to light yellow coloured, crystalline to fine powder
Solubility : Soluble in 1N NaOH
Melting Range : 226 - 236°C
Assay (NT/HPLC) : min. 98.00 % max. 101.00 %

GHS Safety Information

Hazard Statement(s) : H315-H319-H335
Precautionary Statement(s) : P261-P305+P351+P338
Signal Word : Warning
Hazard Pictogram(s)



Risk and Safety Information

R-Phrase(s) : 36/37/38
S-Phrase(s) : 26-36
WGK : 2
RTECS : AU6252200
Storage Temperature(°C) : Store below 30°C
Hazard Symbol(s)



Transport Information

Marine Pollutant : No
ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods