

PGE0011**HiGroNA Powder****Product Number****Packing**

PGE0011 : 50G
PGE0011 : 100G

Product Information

Product Code : PGE0011
Product Name : HiGroNA Powder
Shelf Life : 2 years

Technical Specification

Appearance : White to off-white coloured powder
Solubility : Soluble in 1N NaOH
Melting point : 130 - 133°C
2-Naphthalene acetic acid : <= 4.0%
Assay : min. 95.00 %
Hazard Pictogram(s)

**Risk and Safety Information**

R-Phrase(s) : 22-37/38-41
S-Phrase(s) : 22-26-36
RTECS : QJ 0875000
Storage Temperature(°C) : Store below 30°C
Hazard Symbol(s)

**Transport Information**

ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods