

---

## HiGroIA Powder

**PGE0013**

---

### Product Identifier

HS Code : 2933 99 90  
Storage : Below 30°C  
Shelf life : 2 years

---

### Technical Specification

Appearance : White to yellow to brown crystals or powder  
Solubility : Soluble in 1N NaOH / KOH  
Melting range : 165 - 169°C  
Assay (NaOH Titration/HPLC) : 99.00 - 102.00%

---

### Safety Information

UN No. : Not dangerous goods  
Class : -  
Packing Group : -  
RTECS : Not available  
WGK : 3