

PGE0015

HiMaG Pure

Product Number

Packing

PGE0015 : 500G
PGE0015 : 1KG

Product Information

Product Code : PGE0015
Product Name : HiMaG Pure
Shelf Life : 2 years

Technical Specification

Appearance : White crystals or powder
Solubility : Soluble in water
pH (5 % in water at 25°) : 5.00 - 8.00
Water Insoluble Matter : <= 0.005%
Calcium (Ca) : <= 0.01%
Heavy metal (as Pb) : <= 0.0005%
Iron (Fe) : <= 0.0002%
Manganese (Mn) : <= 0.0005%
Phosphate (PO₄) : <= 0.001%
Sodium (Na) : <= 0.005%
Zinc (Zn) : <= 0.001%
Potassium (K) : <= 0.005%
Copper (Cu) : <= 0.0005%
Chloride (Cl) : <= 0.0005%
Assay (EDTA Titration) : min. 99.50 % max. 100.50 %

Risk and Safety Information

WGK : 1
RTECS : OM4508000
Storage Temperature(°C) : Store below 30°C

Transport Information

Marine Pollutant : No
ADR/RID : Not Dangerous Goods
IMDG : Not Dangerous Goods
IATA : Not Dangerous Goods